

FYD-1 TABLET HARDNESS TESTER HUANGHUA FAITHFUL INSTRUMENT CO.,LTD.

- Tablet hardness testers are instruments for detecting breaking hardness of tablet.

Applicable standards:

Corporate standard (Tablet hardness tester)

Features:

- The tablet is loaded and pressed artificially.
- There is high accuracy press sensor and hardness data is displayed on LED.
- It tests tablet hardness continuously.
- Test data is displayed and locked automatically, so can test be reset and cycled.
- Automation: auto-test, auto-diagnose, auto-alarm.

Specifications:

- Hardness range (2~199.9) N
- Hardness accuracy $\pm 1.5\%$
- Probe procession range 20 mm
- Power 220V/ 50Hz/ 60W or 110V/ 60Hz/ 60W
- Dimension (40*24*16) cm³
- Weight 8 kg



FYD-2 TABLET HARDNESS TESTER HUANGHUA FAITHFUL INSTRUMENT CO.,LTD.

- Tablet hardness testers are instruments for detecting breaking hardness of tablet.

Applicable standards:

Corporate standard (Tablet hardness tester)

Features:

- The tablet is loaded artificially and pressed automatically.
- There is high accuracy press sensor and hardness data is displayed on LED.
- It tests tablet hardness continuously. The biggest number is 100 pieces of a gang.
- Time of testing will be saved by directly measuring the diameter.
- Test data is displayed and locked automatically, so can test be reset and cycled.
- It can deal with and print out such data as: hardness of each tablet, the maximum value, the minimal value, the average value and the number of tablets.
- The parameters can be preset at any time. Either preset or real data is displayed on the LED.
- Automation: auto-test, auto-diagnose, auto-alarm.

Specifications:

- | | | | |
|--------------------------|-------------|----------------------------|------------------------------------|
| ● Hardness range | (2~199.9) N | ● Max number in each group | ≤100 P |
| ● Hardness accuracy | ±1.5 % | ● Power | 220V/ 50Hz/ 60W or 110V/ 60Hz/ 60W |
| ● Probe procession range | (3~40) mm | ● Dimension | (50*30*16) cm ³ |
| ● Test period | ≧ 20 S / P | ● Weight | 12 kg |



FYD-3 TABLET HARDNESS TESTER HUANGHUA FAITHFUL INSTRUMENT CO.,LTD.

- Tablet hardness testers are instruments for detecting breaking hardness of tablet.

Applicable standards:

Corporate standard (Tablet hardness tester)

Features:

- The tablet is loaded artificially and pressed automatically.
- There is high accuracy press sensor and hardness data is displayed on LED.
- It tests tablet hardness continuously. The biggest number is 100 pieces of a gang.
- Time of testing will be saved by directly measuring the diameter.
- Test data is displayed and locked automatically, so can test be reset and cycled.
- It can deal with and memorize such data as: hardness of each tablet, the maximum value, the minimal value, the average value, the number of tablets.
- The parameters can be preset at any time. Either preset or real data is displayed on the LED.
- Automation: auto-test, auto-diagnose, auto-alarm.

Specifications:

- | | | | |
|--------------------------|-------------|----------------------------|------------------------------------|
| ● Hardness range | (2~199.9) N | ● Max number in each group | ≤100 P |
| ● Hardness accuracy | ±1.5% | ● Power | 220V/ 50Hz/ 60W or 110V/ 60Hz/ 60W |
| ● Probe procession range | (3~40) mm | ● Dimension | (40*24*16) cm ³ |
| ● Test period | ≧ 20 S / P | ● Weight | 8 kg |

